

DRAMINŃSKI *Animal profi2*

UNIVERSAL VETERINARY ULTRASOUND SCANNER FOR LARGE AND SMALL ANIMALS

- possible to connect different types of mechanical probes
- exchangeable effective battery
- solid frame body
- sun protection shield



Call:
+48 89 675 26 00

Write:
info@draminski.com

DRAMIŃSKI *Animal profi2*

**UNIVERSAL VETERINARY ULTRASOUND
SCANNER FOR LARGE AND SMALL ANIMALS**



- ▶ It is possible to connect different types of mechanical probes to the scanner for rectal and abdominal tests,
- ▶ It has an exchangeable effective battery which enables 4.5 hours of continuous work in the field,
- ▶ It has a solid frame body which guarantees safety during operation, it is easy to clean and disinfect after each test,
- ▶ It is possible to attach a sun protection shield which makes working in direct sunlight easy.



DESIGNED TO WORK IN THE FIELD

Standard package includes: an ultrasound scanner with a selected mechanical probe, 2 external batteries, a charger, USB adapter, neck and hip strips, ultrasound gel (250 mm), transport case, user manual.

Technical specification:

- Size 23.0 x 17.0 x 6.0 cm
- Weight 1530 g
- 7" LED-backlit LCD display

Image:


- B, B+B, B+M,
- post-processing – digital image processing,
- zoom 60 – 200%,
- freeze,
- regulation: gain, gamma, frequency.

Probes (mechanical, sector):


- rectal 3 – 7 MHz, scan angle 180°, penetration depth up to 12 cm,
- abdominal 3 – 7 MHz, scan angle 90°, penetration depth up to 20 cm,
- rectal-abdominal "Down Fire" 3 – 7 MHz, scan angle 90°, penetration depth up to 20 cm,
- rectal with a holder 3 – 7 MHz, scan angle 180°, penetration depth up to 12 cm.

Functions:

- save images and cine loop with description and date to the ultrasound scanner's disc,
- transmission of images and cine loop to external data carrier,
- measuring (distance, surface, volume, fetal biometry charts),
- regulation of scanning range,
- full screen.



Sow – pregnancy



Cow – ovary

DRAMIŃSKI S.A.

Wiktora Steffena 21, 11-036 Sząbruk, Poland
tel. +48 89 675 26 00, ultrasound@draminski.com

**WE CONSTRUCT
WE PRODUCE
WE SELL
WE SERVICE**