



SEE MORE

DRAMIŃSKI S.A.

Wiktora Steffena 21, 11-036 Sząbruk, Poland

+48 89 675 26 00 @ export@draminski.com

Follow us:    

www.draminski.com

FOR THE BEST
LIVESTOCK BREEDERS

PERFECT FOR FARM WORK

ESTROUS DETECTORS


The estrous detector helps breeders to determine the exact moment of the estrous cycle when the ovulation takes place (also in the case of silent heat).

Main features:

- large display with energy-saving LED backlighting
- safe for animals and humans
- does not require periodic calibration
- memory for 250 animals and 200 thousand measurements with date and time
- 2-year warranty



The estrous detector is available in version for:
cows and mares | sows | sheep and goats



Examples of typical changes in electrical resistance of vaginal mucus registered by the DRAMINSKI estrous detector



PREGNANCY DETECTORS

Ultrasound Pregnancy Tester is invaluable help in modern breeding. The detector provides maximum help for breeders, advisers and veterinarians in making an accurate diagnosis that confirms or excludes pregnancy.

Main features:


- early diagnosis
- easy to use
- high efficiency rate (almost 100%)
- easy disinfection

The estrous detector is available in version for:
sows | sheep and goats



Animal species	The days of examination from-to	The optimal day of examination
pig	18 - 75	25 - 30
sheep	50 - 100	70
goat	40 - 100	60 - 70

Demonstration of correct application of the probe and optimal test time



Examples of typical changes in the electrical resistance of milk

MASTITIS DETECTORS

The device to detect subclinical mastitis in its initial and invisible stage.

Main features:

- large display with energy-saving LED backlighting
- safe for animals and humans
- does not require periodic calibration
- memory for 250 animals and 200 thousand measurements with date and time
- USB socket
- 2-year warranty

FOR THE MOST RELIABLE RESULTS THE FIRST JETS OF MILK SHOULD BE USED DIRECTLY IN THE MEASUREMENT CUP

ULTRASOUND SCANNERS



ANIMAL PROFI 2

Portable ultrasound scanner to control animal's reproduction process.

Main features:

- post-processing
- small dimensions
- long battery life
- perfect image even in direct sunlight
- data transfer to PC



SONOFARM MINI

Portable ultrasound system for rapid pregnancy detection.

Main features:

- perfect for field work
- interchangeable probes
- rapid diagnosis
- easy to use
- dedicated to breeders, inseminators and veterinarians



- 1 Sector rectal 5,0 MHz, 180°, on extender
- 2 Sector rectal 5,0 MHz, 180°, on extender
- 3 Abdominal down-fire 5,0 MHz (sector) ø31 mm
- 4 Sector rectal 5.0 MHz, 180°
- 5 Abdominal (sector) 5,0 MHz

