




airScan
Boost your efficiency

Dramiński **airScan** is a **lightweight, wireless ultrasound scanner** for farm animals diagnostics.

Thanks to the latest technology, it provides **superior image quality**.

The **Colour Doppler** mode will increase your efficiency.

Available in two versions - with a **Rectal Linear** or **Rectal Convex** probe.



Application

	airScan Rectal LINEAR	airScan Rectal CONVEX
Bovine reproduction	+	+
Equine reproduction	+	+
Small ruminant reproduction		+
CL blood flow	+	
Lungs	+	+
Body condition	+	+
Udder and teat	+	
Liver		+
Gastrointestinal tract	+	+
Tendons	+	

Key features

- Includes **Colour Doppler** mode
- Weighs **only 0,6 kg**
- Short boot up time
- Quick and easy probe connection
- Up to **4,5 h of working time** on one charge (B mode)
- Compatible with Android and iOS system
- **No subscription required**

Imaging modes

- B mode
- B+B mode
- Colour Doppler mode
- Power Doppler mode

Measurements

- Grid
- Distance
- **Age of foetus**

Standard set up contains

- **airScan** Rectal Linear or Rectal Convex
- Quick mount
- Waist belt
- Charger
- Transport case

Hardware requirements

- **iOS**
iOS 11 or above
iPhone 7 or above
iPad 3 or above
- **Android**
Android 7.0 or above
CPU: Snapdragon 650 or above
RAM: 4GB or above

airScan



Dramiński airScan App
Connect the airScan probe to your device.





SEE MORE

DRAMIŃSKI S.A.

Wiktora Steffena 21, 11-036 Sząbruk, Poland

 +48 89 675 26 00  ultrasound@draminski.com

Follow us:    

www.draminski.com