



iScan 2 MULTI

comprehensive diagnostics



Ultrasound scanner **iScan 2 MULTI**

iScan 2 MULTI is a multi-task field ultrasound scanner with interchangeable probes for comprehensive diagnostics of small and large animals.

NEW IMAGE QUALITY

LuciD™ is a technology reducing the noise and increasing the contrast, sharpness and tissue differentiation in the image. It makes the image clearer and easier to interpret.

LARGE SCREEN

7" LCD screen with wide viewing angles ensures comfortable work. We know that you often face difficult conditions. That is why we rugged the screen making it less prone to damage.

START EXAMINATION IN A FEW SECONDS

Adjust the scanner to your needs - select the most used functions and save them in the quick access menu (Q). These one-button access will make your work faster and more friendly.

CLEANING AND DISINFECTION

Flat screen, probe socket protection slide and waterproof keys make cleaning the ultrasound scanner easier and more comfortable.

STURDY AS ALWAYS, LIGHTER THAN YOU EXPECTED


Rugged casing makes the **iScan 2 MULTI** extra durable. The active ventilation system implemented in the casing ensures good operation parameters regardless of heat and humidity.

HIGH PERFORMANCE, NON-STOP

A powerful battery provides up to 7 hours of non-stop operation. In addition, you can turn off the screen to save power without shutting the scanner down. And the included suspenders will not let you get tired for the entire day.

dimensions: 256 x 165 x 69 mm (10.1 x 6.5 x 2.7 in)

weight: 2.06 kg (4.54 lbs)



probe socket protection slide

operating 7 h on one battery



Set contains

- **iScan 2 MULTI** ultrasound scanner
- Ultrasound probe of your choice*
- External battery
- Battery charger
- USB adapter
- Carrying straps
- Transport case
- User guide (USB)

* Select the probe prior to purchase

Wide range of available probes



Lots of probes to choose from: 1 Ovum Pick-Up System 2 Microconvex for OPU 3 Convex 3.5 MHz 4 Convex 5.0 MHz 5 Linear rectal 6 Convex rectal 7 Backfat.

Accessories



OLED goggles




Sun protection shield



Extension for linear rectal probe




WiFi antenna




Android application

DRAMIŃSKI
CARE AND SHARE **app**


Diagnostic images




Endometritis



Inflammation of the umbilicus



Cyst



Pregnancy - day 28



SEE MORE

DRAMIŃSKI S.A.

Wiktora Steffena 21, 11-036 Sząbruk, Poland



+48 89 675 26 11



ultrasound@draminski.com

www.draminski.com