

Double Up Your Efficiency
XCUBE i9



ALPION
MEDICAL SYSTEMS

Proven **premium solutions** – now tailored for veterinary medicine

+48 89 675 26 00

@ sales@draminski.com

Double Up Your Efficiency

XCUBE i9

The **Alpinion XCUBE i9** offers uncompromising performance and portability. With its excellent image quality and advanced features, it offers veterinary specialists the highest level of diagnostics.

What sets us apart

Discover selected highlights of the new veterinary-dedicated ultrasound line:

- **Configuration tailored to your needs**
choose the setup that matches your caseload and workflow
- **Cardiology module upgrade at any time**
the system grows with your practice, whenever you're ready
- **Reinforced linear probe**
improved durability for intensive daily use
- **2 active probe ports**
fast probe switching without re-plugging or downtime



Optional features

- **X+MicroView** | sensitive microvascular Doppler for imaging very small-vessel flow
- **X+Compare** | side-by-side comparison of images from different exams of the same patient
- **X+Assistant** | exam protocols supporting repeatability and standardisation
- **Elastography** | strain elastography for linear probes
- **Warranty up to 5 years** (optional) | standard warranty: 2 years

Probes

- **SC1-7H** | Wideband convex probe (1-7 MHz), 192 elements, X+ Crystal Signature technology
- **SC2-11H** | Wideband micro-convex probe (2-11 MHz), 192 elements, X+ Crystal Signature technology
- **SL3-19H** | Single-crystal wideband linear probe (3-19 MHz), 40 mm footprint, 192 elements
- **L10-25H** | High-frequency wideband linear probe (10-25 MHz), 25 mm footprint, 192 elements
- **MP1-5X** | Multi-layer phased array probe (1-5 MHz), X+ Crystal Signature technology
- **SP3-8T** | Single-crystal phased array probe (3-8 MHz)