




v.198712125EN

**DRAMIŃSKI**  
TECHNOLOGY

www.draminski.com



SEE MORE

**DRAMIŃSKI S.A.**

Wiktora Steffena 21, 11-036 Sząbruk, Poland

+48 89 675 26 00 @ ultrasound@draminski.com

Follow us:

[www.draminski.com](http://www.draminski.com)

**iLens**  
**WIRELESS GOGGLES**  
FOR PRECISE ULTRASONIC DIAGNOSTICS IN ANIMALS

# iLens

## WIRELESS GOGGLES



Elevate your ultrasound experience to unprecedented heights with iLens. Witness the revolutionary impact of AR technology transforming ultrasound examinations.

**Introducing iLens:** cutting-edge wireless goggles meticulously crafted for veterinarians. Powered by advanced augmented reality (AR) technology, these goggles offer unparalleled freedom and comfort, revolutionizing the ultrasound examinations. Experience a new pinnacle of efficiency with iLens.

### Benefits:

**Wireless, seamless image transmission** with full refresh rate speed.

**Always clear image** - as well as on those sunny days outdoors.

Experience **the utmost safety** when working with animals, courtesy of state-of-the-art AR technology.

Enjoy an **expansive field of view** and complete control over your surroundings and movements during examinations.

Ability to work with **prescription glasses**.

Full-day autonomy. **Over 12 hours** on a single battery charge.

Ability to **charge in the car** and work with a power bank\*.

Sturdy headpiece and smooth circumference adjustment for **secure wearing on the head**.

Perfectly positioned center of gravity ensures smooth movement and **reduces strain on the neck and head**.

No cables guarantee **greater freedom**.

**Easy pairing and connection** with the ultrasound machine.

**Dedicated carrying case** with space for a charger, USB-C cable, sunshade, and accessories.

### PRECISION, MOBILITY SAFETY



### Specification:

- Display: transparent
- Resolution: 854 x 480
- Visor: transparent, dimmable
- Connectivity: Wi-Fi 802.11ac
- Ports: USB-C
- Battery: 6000 mAh
- Weight: 465 g

### Components:

- 1 iLens goggles
- 2 Carrying case
- 3 Bright screen visor
- 4 Dark screen visor
- 5 USB-C cable
- 6 USB charger
- 7 Power bank\*



**iLens goggles are compatible with selected ultrasound machines, please inquire about details from the distributor.**



\*Power bank is not included in the set. It is available as an additional option.